







在正式使用设备前 ...



称量



技术参数

200 kg	2 kg	$< 150 \text{ kg: } 100 \text{ g}$ $\geq 150 \text{ kg: } 200 \text{ g}$	 	6 x AA	+5 °C ... +40 °C	$a = 321 \text{ mm}$ $b = 356 \text{ mm}$ $c = 60 - 66 \text{ mm}$	4.0 kg

校准：按照等级 III 经过医学校准，符合 93/42/EEC 指令 I 级医疗产品要求

保持
带有呼叫标记的三角图标表示，显示冻结时涉及的是无需校准的辅助装置。

保修说明

对于由设备本身或制造引起的质量问题，我们将为您提供自供货起的两年保修期。所有可活动部件如电池、电缆、电源设备、蓄电池等除外。您只需出示购买发票，我们便将为您免费解决一切属于保修范围内的质量问题。其他的问题不在我们的责任范围之内。如设备所在地与客户所在地不一致，则因此产生的来回运输费用由客户承担。

对于由运输引起的设备损害，我们只在下列情况下承担保修责任：运输设备时使用了全部的原包装且设备按照原包装状态安全固定好。因此，请您保管好所有的包装件。如设备由非赛康授权人员打开，我们将不承担保修责任。如国外的客户有保修需要，请直接与当地的销售商联系。

地址

赛康医用度量系统（杭州）有限公司
萧山经济技术开发区桥南区鸿达路 201 号
浙江杭州（邮政编码 311231）
中国

电话: +86 571 82 86 96 89
传真: +86 571 82 86 96 87
邮箱: info.nlv@seca.com

Medical Measuring Systems and Scales since 1840

seca gmbh & co. kg
Hammer Steindamm 3-25
22089 Hamburg · Germany
Telephone +49 40 20 00 00 0
Fax +49 40 20 00 00 50
info@seca.com

seca operates worldwide with headquarters
in Germany and branches in:

- seca france
- seca united kingdom
- seca north america
- seca schweiz
- seca zhong guo
- seca nihon
- seca mexico
- seca austria
- seca polska
- seca middle east
- seca brasil
- seca suomi
- seca américa latina
- seca asia pacific

and with exclusive partners in
more than 110 countries.

All contact data under www.seca.com

